



Key4careTM

Key4care solves the problem of property access for community or health-care workers and emergency staff in an effective and cost efficient manner. It promotes independent living as part of a planned care package or as an integral part of a social alarm program.

- ◆ Improved response time to emergency calls
- ◆ Reduced costs of key management
- ◆ Easily transferred to new locations
- ◆ Suitable for use with deliveries of supplies



When you're dealing with health-care issues in people's homes, gaining access can be a significant problem. Finding a keyholder, diverting to pick-up a key from a central location or using a fixed access code, are just not quality solutions.

The Key4care digital keybox is an ideal alternative:

- ◆ **It is a discrete and tough little box bolted to a wall**
- ◆ **It holds a key and requires an access code for entry**
- ◆ **Access codes are provided "live" when entry is needed**
- ◆ **Each box is unique and each code is valid for only a short period of time**

Approx dimensions:
105mm x 60mm x 140mm
Powered by 4 x AA (LR6)
batteries

It's easily installed - 3 holes is all it takes. There is no wiring to mains power or telephone lines. It can also be easily moved to a new location when required.

You are logged into 'UK', where the current time is 24/02/2005 10:59:32 (UTC+1)

Key Generation

Date/Time:	02 / May / 2005 10 : 59 : 02	Auto: <input checked="" type="checkbox"/>	?
Unit serial:	075B	or select: <input type="text"/>	?
	Lookup address		
House name/number:	111		?
Street:	Lime Street		
Postcode:	TA23 9AS		
	Find address		
Mode of use:	<input checked="" type="radio"/> Normal <input type="radio"/> Demo		?
CGS mode:	Secure Access		
Access restrictions:			?
Caller ID:	7		?
	Validate caller		
Caller name:	Wally Smith "favourite curry"		
Security question:	Please enter 2nd letter: <input type="text" value="h"/> Please enter 7th letter: <input type="text" value="k"/>		?
Get key code:	Get access code	Restart	
Authorised key code:	3652		

Managed on-line

- ◆ All details on the estate of keyboxes can be managed online (names, addresses, serial numbers, authorised users)
- ◆ Access codes are generated "live" and details recorded for traceability
- ◆ As employees leave, they can simply be removed from the list of authorised users